Innovations in orthodontics


Gebrauchsanweisung
User Manual

PowerScrew 8, 12, 16

Dehnschraube zur Gaumennahterweiterung
Expansion screw for maxillary expansion

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
Notizen

A large grid of small dots for taking notes, consisting of 20 columns and 25 rows.

English

The components of the PowerScrew

Technical information

| Screw bodies (size) | Max. expansion path | Ring gauge |
|---------------------|---------------------|------------|
| 8 mm | 6 mm | 3 rings |
| 12 mm | 12 mm | 5 rings |
| 16 mm | 18 mm | 7 rings |

Expansion per turned surface

| Surfaces | expansion path |
|----------|-------------------|
| 6 | 1,0 mm (= 1 turn) |
| 1 | 0,17 mm |

Dear doctor, we are honored that you have chosen this innovative new product from TigerDental.

To ensure the safe and effective use of this product, you must carefully read and adhere to this user manual.

The user manual covers all necessary safety instructions but cannot describe all possible circumstances. For this reason, you may find further user notices online at tigerdental.com under the "Know-how" link. Our customer service is available for any additional questions.

The continual development of our products indicates our company's commitment to innovation. Accordingly, after several clinical uses, we recommend you carefully review the accompanying user manual (this can be found online under the "Downloads" link at tigerdental.com).

1. Manufacturer

TigerDental, Belruptstraße 59, 6900 Bregenz, Austria

2. Product description

The PowerScrew is a double telescope screw for use in palatal expansion devices. It is available in three sizes (8, 12 and 16 mm installation size).

Activation is simple, as the 8 cm long flat spanner can be quickly and safely used by the patient on the hexagonal screw.

Narrow palates need a small expansion screw. The PowerScrew therefore grows with it. By simply changing the hexagonal screw for the next biggest, the expander can continue to be used without making a new one (see images page).

The PowerScrew and its retention arms are extremely stable and can withstand the greatest forces and torques.

One turn (6 surfaces) equals 1 mm expansion (0.17 mm per surface); every second surface is marked with increasing points..

3. Range of application

Indications: Expansion of skeletally narrow upper jaws

Contraindications: Age-dependent, ossified palatal expansion

4. Use of the Expander and activation of the PowerScrew

As soon as the device is used, medical activation of at least 0.5 mm (3 surfaces) can be carried out. That guarantees initial activation and shows the patient the correct use.

Please give your patient the patient information leaflet with the turning protocol for daily recording. Enter in this turning protocol the turns made by you for activation. We also recommend allowing the patient to make a few turns forwards and back for him/her self. This ensures that the patient has understood the turning process.

Please carry out a withdrawal test at every checkup. Use the Weingart forceps to pull briefly on the retention arms to check the seat and firmness of the apparatus.

5. Expansion, warning thread and turning stop

Each ring gauge corresponds to 2 mm expansion. When maximum expansion of the PowerScrew is reached, the warning thread (green) appears on the threaded tube. The screw should not be turned any further.

6. Securing the PowerScrew after turning stop

After expansion is concluded, light-cured resin (Heliosit® or Triad® Gel) can be applied between the hexagonal screw and the retention arms to fix the hexagonal screw to the thread rod. Further expansion can be carried out after that if required.

7. Instructions for single use products

The PowerScrew is ready to use in its packaging but it is not sterile and should only be used (installed) once. Repreparation of a single use PowerScrew (recycling) or further use on the patient is not allowed.

8. Safety instructions

The above guidance for use and implementation includes tips and recommendations, but this does not release the doctor from personal responsibility!

According to the law on medicinal products (MPG), the use of the PowerScrew has been approved exclusively for trained orthodontists and dentists. The user is expressly instructed to adhere to its use in accordance with the terms of the law on medicinal products and any other applicable regulations and laws. The user himself/herself is responsible for compliance.

Patients with an allergy to nickel or chrome should not use this product. If there is any irritation the advice of a doctor must be sought immediately for interruption of treatment.

The PowerScrew must be removed before magnet resonance tomography (MRI) is carried out. As soon as treatment is concluded, the PowerScrew can be implemented again.

9. Quality instructions

TigerDental assures the user of the faultless quality of the product. The contents of these instructions for use are based on self-experience. The user is responsible for the processing of the product itself. In default of exertion of influence by TigerDental on processing by the user, there is no liability for faulty results.

Although the PowerScrew is a mature product, improvements and changes may be carried out which may not be recorded in the above instructions for use. Nor do we accept any liability for printing errors.

10. Contents

- 1× PowerScrew with/without retention arms
- 1× flat spanner

11. Description of the used symbols

A description of the used symbols of ÖNORM EN ISO 15223-1 can be found at: www.tigerdental.com/downloads

By changing the hexagonal screw up to 18 mm expansion possible with one device

The PowerScrew 12 and 16 with retention arms

Examples of palate expansion apparatus

Cap splint device

Haas device

Hybrid device with mini anchor screws

MICRO2 (purely skeletal, 12–16 years)

MICRO4 (purely skeletal, from 17 years of age)

Simple activation with flat spanner by the patient himself/herself

Example of treatment

Eight-year-old patient before palate expansion

Cap splint device with PowerScrew12 11 mm expansion, typical diastema

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